


## Measurement of pathway-specific signaling kinetics in real-time with bioSensAll®

Comprehension of the temporal aspects (i.e., kinetics and dynamics) of receptor-mediated signaling is integral to fully appreciating the mechanisms by which extracellular signals impact cellular behavior. For ligand-receptor interactions, kinetics can be generally described at two levels: i) ligand binding to the target receptor, and ii) receptor downstream signaling. Ligand binding kinetics are defined in terms of  $k_{on}$  and  $k_{off}$  values (ligand association and dissociation rates, respectively) determined via radioligand binding assays. Signaling kinetics are defined in temporal terms and may range from fractions of a second to many hours.

The effect of kinetics on a ligand's efficacy and biological action(s) has long been recognized but all too often overlooked [1]. Consideration of kinetics is all the more important when dealing with G-protein coupled receptors (GPCRs) given their well-appreciated capacity to activate multiple signaling pathways exhibiting distinct kinetic properties. Further, it is now evident that different ligands of a given GPCR can activate distinct subsets of

receptor-coupled signaling cascades, thus leading to compound-specific signaling signatures (a.k.a., biased agonism). Recent observations have suggested that both a ligand's signaling signature and biased qualities may be qualitatively and quantitatively modulated over time [2]. Single point readouts of signaling conducted at equilibrium may therefore provide an incomplete portrait of a ligand's character and fail to capture important differentiating elements between compounds. As a result, the ability to measure signaling kinetics in real-time over short (milliseconds) and extended (hours) timescales offers invaluable insight into a compound's pharmacological properties.

In this application note, we demonstrate how a single, multiple-compound addition assay with **bioSensAll®** biosensors was used to obtain real-time kinetic measurements of G-protein (GoA) activation and  $\beta$ -arrestin 2 membrane recruitment for the human mu-type opioid receptor (hMOR) agonist and antagonist DAMGO and naloxone, respectively.



**Results and Conclusion** – DAMGO-induced GoA activation and plasma membrane recruitment of  $\beta$ -arrestin 2 displayed comparable kinetics post-agonist addition. However, whereas maximal GoA activation was sustained for at least 30 min, the presence of  $\beta$ -arrestin 2 at the plasma membrane declined in a time-dependent manner. Pathway-specific differences in sensitivity to the antagonist activity of naloxone were also revealed. Indeed, naloxone displayed faster kinetics in antagonizing DAMGO-

dependent  $\beta$ -arrestin 2 membrane localization (especially visible with naloxone IC<sub>50</sub>). **bioSensAll®** technology allows for the real-time assessment of ligand signaling kinetics and dynamics across multiple GPCR-downstream signaling pathways and on timescales ranging from milliseconds to hours in a simple and homogenous assay format. **bioSensAll®** thus provides additional information about ligand pharmacology that cannot be obtained via equilibrium or endpoint assays.

### References

[1] Mackay, D. (1977) A critical survey of receptor theories of drug action. *Kinet. Drug Action*. 47:225–321. [2] Klein Herenbrink, C., Sykes, D.A., Donthamsetti, P. et al. (2016) The role of kinetic context in apparent biased agonism at GPCRs. *Nat. Commun.* 7. Article number: 10844.

more information: [collabdtna@domaintherapeutics.com](mailto:collabdtna@domaintherapeutics.com)

more on our website: [domaintherapeutics.ca](http://domaintherapeutics.ca)