## HIGH-THROUGHPUT APPLICATIONS FOR SCREENING GPCR ACTIVATION USING DIOSensAll® THERAPEUTICS

Instrumentation



Frank Cao<sup>1</sup>, Claire Normand<sup>1</sup>, Arturo Mancini<sup>1</sup>, Samya Aouad<sup>1</sup>, Herthana Kandasamy<sup>1</sup>, Robert Cerone<sup>1</sup>, Meriem Semache<sup>1</sup>, Melissa R. Miller<sup>2</sup>, Shunjie Guan<sup>2</sup>, Dong Yan<sup>2</sup>, Craig L. Hyde<sup>2</sup>, Jean-Philippe Fortin<sup>2</sup>, Eric B. Fauman<sup>2</sup>, Stephan Schann<sup>1</sup>, and Laurent Sabbagh<sup>1</sup>

<sup>1</sup>Domain Therapeutics, Montreal, Quebec, Canada, and Strasbourg, France and <sup>2</sup> Worldwide Research, Development and Medical, Pfizer Global R&D, Cambridge, MA, USA

**Microplate Reader** 

## Abstract

G protein-coupled receptors (GPCRs) represent the largest class of membrane-bound proteins which regulate physiological processes by engaging diverse signaling pathways. Although dysregulation of GPCR activity has been linked to pathological consequences, the involvement of specific signaling pathways toward distinct clinical repercussions remain unresolved. Thus, interrogating the complex signaling behaviors of GPCRs is not only important for understanding GPCR biology, but also essential for successfully targeting these receptors as therapeutic avenues. Here, we describe the application of our bioSens-All® technology within the context of a high throughput screening platform as well as a mechanical stimulation approach to characterize GPCR signaling.

bioSens-All<sup>®</sup> consists of a panel of 16 selective enhanced bystander bioluminescence resonance energy transfer (ebBRET)based biosensors that monitor the activation of heterotrimeric G proteins and β-arrestins upon GPCR activation. To achieve highthroughput volumes, bioSens-All<sup>®</sup> assays were miniaturized and adapted to a 384-well format, where HEK293 cells were cotransfected with receptors and biosensors. Assays were performed using a laboratory automation workflow software to coordinate the interaction between various instruments via a robotic microplate mover. Through this automated platform, we assessed signaling profiles of ~12,000 clinically relevant variants across 60 GPCRs. Analysis was performed through an automated tool to generate and compile over 120,000 concentration-response curves. Signaling profiles for a subset of variants revealed gain- or loss-of-function properties across different signaling pathways thus also allowing us to identify genetic alterations that impart biased signaling profiles for each profiled GPCR.

Using the bioSens-All<sup>®</sup> technology in a separate application, we developed a novel *in vitro* assay detecting adhesion GPCR ADGRE5 activation following mechanical stimulation. In this assay, cells co-transfected with receptors and biosensors were seeded in 96-well plates and subjected to vigorous orbital shaking prior to BRET measurements. We demonstrated that mechanostimulation of full-length human ADGRE5-transfected cells resulted in a rapid time-dependent recruitment of β-arr2 to the plasma membrane, while such mechanosensitivity was not observed with a control receptor (mu opioid receptor, MOR) despite MOR's ability to engage β-arr2 following ligand stimulation. Finally, the assay demonstrated a Z' value of 0.73, thus highlighting its compatibility for high-throughput screening of novel therapeutic agents acting to modulate physiologically relevant modes of ADGRE5 activation.



**Ligand Dispenser** 







Altogether, these applications of our bioSens-All<sup>®</sup> technology allows interrogation of diverse activation mechanisms and complex signaling behaviors of GPCRs. Ultimately, such approaches can be applied to other GPCRs to facilitate the validation of **BRET Substrate Injector** novel therapeutic targets and discovery of innovative therapeutic agents.

βarr **β**-arrestin engagement PROXIMAL

Fig. 1. BRET-based biosensors for G protein and Barr activation. (a) The bioSens-All® platform consists of 16 biosensors for monitoring activation of specific Ga proteins and Barrs downstream of GPCRs. These biosensors measure receptor proximal events engaged upon receptor activation, namely G protein activation and  $\beta$  arr recruitment. In addition to proximal biosensors, the technology includes sensors detecting distal effector activity and second messenger generation. (b) Upon receptor activation, Rluclltagged effector proteins (Effector-RlucII) translocate towards and interact with active Gα subunits from each G protein family, leading to increased BRET. (c) Upon receptor activation, RlucII-tagged βarrs (βarr-RlucII) translocate to the plasma membrane, thus increasing BRET with membrane-anchored rGFP.

## **Mechano-Stimulation Assay**



## **Automated Screening Platform**



**Automated Incubator** 

Fig. 3. Signaling profiles of variant receptors. HEK293 cells were transfected with plasmids coding for a specific GPCR and one of 16 biosensors. Cells were treated with GPCR-specific ligand and response recorded 30 mins later. (a) Variant displays WT-like profile on all pathways selected. (b) Variant displays complete loss of function profile on all pathways selected. (c) variant displays loss of function on G15 and WT-like profile on Gi1, Gi2, and βarr2, thus yielding an overall biased profile.

![](_page_0_Figure_22.jpeg)

Fig. 6. Measuring mechano-stimulation of hADGRE5 through β-arr Fig. 7. Z' robustness of mechano-stimulation BRET assay. plasma membrane recruitment sensor. (a and b) HEK293 cells co- HEK293 cells co-transfected with β-arr plasma membrane transfected with biosensor and hADGRE5 or uncleavable GPS point recruitment sensor and hADGRE5 or MOR were subjected or not

![](_page_0_Picture_25.jpeg)

Fig. 4. Activity of 265 variants of ADRB1. (a) Location of variants on the serpentine map of ADRB1. (b) In vitro pharmacology of G15 area under the curve (AUC) for all profiled ADRB1 variants. (c) G15 signaling of four representative variants.