

# DEVELOPMENT OF A BRET BIOSENSOR-BASED SCREENING STRATEGY FOR IDENTIFYING DRIVERS OF HALLUCINOGENIC ACTIVITY OF 5HT2AR TARGETING COMPOUNDS


Guilhem DUGAST, Arturo MANCINI, Claire NORMAND and Laurent SABBAGH  
Domain Therapeutics NA Inc., Montreal, Québec, Canada

## INTRODUCTION

Promising clinical trial outcomes involving across various indications have stimulated a resurgence of scientific interest in compounds targeting the 5-HT<sub>2A</sub> receptor (5-HT<sub>2A</sub>R). However, many of these molecules act as hallucinogens, hampering their widespread therapeutic use. One goal in this field of research is the development of compounds that retain their therapeutic properties but lack hallucinogenic activity. Yet, there is no clear consensus on the key factors driving the hallucinogenic actions of 5-HT<sub>2A</sub>R agonists. In this work, we describe the application of novel enhanced bystander (eb)BRET-based biosensors towards the characterization and differentiation of hallucinogenic and non-hallucinogenic neuroplastic 5-HT<sub>2A</sub>R targeting molecules. These biosensors allow for the real-time and broad multiparametric monitoring of the signal transduction pathways engaged upon activation of native G Protein-Coupled Receptors (GPCRs). Using these tools, we can evaluate five core pharmacological factors that may help in distinguishing hallucinogenic from non-hallucinogenic neuroplastic 5-HT<sub>2A</sub>R targeting molecules: pathway-specific compound efficacy / signaling bias at 5-HT<sub>2A</sub>R, compound polypharmacology, receptor activation kinetics, receptor internalization, and compound activity at subcellular compartments. We propose a biosensor-based and AI-supported screening pipeline in which various hallucinogenic and non-hallucinogenic 5-HT<sub>2A</sub>R agonists are first profiled on these five pharmacological parameters and the resulting data subsequently used to train machine-learning algorithms to cluster compounds based on their hallucinogenic activity. Novel 5-HT<sub>2A</sub>R targeting compounds can then be funneled through this pipeline to help isolate non-hallucinogenic compounds, and whose activity could then be validated *in vivo*. This approach, used in an iterative manner, may help inform the early identification and rational design of safer, non-hallucinogenic compounds for the use in various neuropsychiatric indications.

## MATERIALS AND METHODS

### ebBRET-based biosensor platform **BioSensAll**<sup>®</sup>




**Figure 1:** Assay principle underlying our ebBRET-based biosensors

Avet, Mancini et al. eLife. 2021.

## RESULTS & CONCEPTUAL FRAMEWORK


### Compound signaling bias & efficacy



**Figure 2:** Multiparametric signaling profiling allowing for a qualitative and quantitative differentiation of compound signaling capacities. Shown here is the signaling profile of the 5-HT<sub>2A</sub>R receptor in response to 5-HT (serotonin).

Comparison of compound activity on various signaling pathways enabling the detection of inter-pathway differences in compound efficacy and functional selectivity. Differences in efficacy and functional selectivity may drive hallucinogenic activity of 5-HT<sub>2A</sub>R targeting agents (Wallach et al. DOI 10.1101/2023.07.29.551106. 2023; Karaki et al. Mol Cell Proteomics. 2014).


### Compound polypharmacology




**Figure 3:** Psychoactive substances often display activity at multiple GPCRs (i.e., polypharmacology). The ability to detect such promiscuity may allow to identify polypharmacological patterns that distinguish hallucinogenic from non hallucinogenic 5-HT<sub>2A</sub>R compounds. Shown here is the activity of various ligands on different receptors using our ebBRET biosensor assays. Avet, Mancini et al. eLife. 2021.

### Receptor activation kinetics

#### Gq activation of h5-HT<sub>2A</sub> receptor by 5-HT (Serotonin)




**Figure 4:** Receptor binding kinetics of different compounds may influence their therapeutic properties and their hallucinogenic potential. LSD's slow binding kinetics translates into a shift in potency. This may explain LSD's long duration of action and hallucinogenic activity *in vivo*.



**Figure 5:** 5HT<sub>2A</sub>R binders can display unique activation kinetics on each G protein and Barrestin pathway. Ligand-specific differences in pathway activation kinetics may translate into differences in hallucinogenic and therapeutic potential.


### Compound activity at subcellular compartment



**Figure 7:** Detection of compound-induced receptor activity at intra-cellular compartments; it has been proposed that psychedelics promote neuroplasticity through the activation of intracellular (Golgi) 5-HT<sub>2A</sub>R (Vargas et al. Science. 2023).

### Receptor internalization

**Figure 6:** 5HT<sub>2A</sub>R internalization is another parameter that can serve to differentiate various binders. Shown is our GPCR internalization assay measuring a time-dependent reduction of cells-surface GPCR levels.



## Non-hallucinogenic

## Hallucinogenic

## CONCLUSIONS

- ebBRET-based biosensors can be applied towards advanced multiparametric pharmacological characterization of 5-HT<sub>2A</sub>R targeting agents.
- Measurement of 5 key pharmacological parameters, coupled with machine learning algorithms, can help quickly distinguish hallucinogenic and non-hallucinogenic 5-HT<sub>2A</sub>R binders.
- The proposed biosensor-based and AI-supported screening pipeline can be used to identify and generate ligands with the desired psychoactive properties.